

Exergen's New Keyless Self-Locking Wall Mount with Hidden Release



134305
Keyless Locking
Wall Mount Instrument
Instrument Holder



134306
Keyless Locking Wall Mount.
Holder and Cap Holder
Holder plus integrated cap
dispenser shown with caps



Security Cable - Field Attachment

8 ft. Coiled Rubber #124309

8 ft. Coiled Rubber - Latex Free Order #124311

6 ft. vinyl-covered steel, Order #134302

8 ft. vinyl-covered steel, Order #134303

All models of the Keyless Locking Wall Mount illustrated above can be used with or without cables. Security cables are sold separately.

Suggestions for Selection:

- If instruments are to be secured at the nurses station, consider the keyless locking wallmount 134305 or 134306 without a cable.
- If instruments are to be mounted one per bed, or in a dedicated area (triage for example), consider 134201 with the coiled cable 124309 or 124311.
- If instruments are to be attached to an existing rolling stand, or a vital signs monitor, consider the vinyl covered steel cable 134302